


Ultherapy

SKIN TIGHTENING SYSTEM



W o r k T h e o r y

- Ultherapy is Non Surgical Face Lift Treatment- Ultherapy (High Intensity Focused Ultrasound) Technology , the aging effects will show on the skin like the wrinkles and the eye-bags. It is the signals to alert you your collagen and elastin are losing from your body. The collagen is the key structural component to support your skin.



- Ultherapy (High-Intensity Focused Ultrasound) which is the new technical for face lifting and tightening.
- Ultherapy is delivered into the deep skin and cause thermal coagulation to SMAS (Superficial Muscular Aponeurotic System). It can provide the skin tightening and rejuvenation happens gradually with the time.
- Ultherapy technology device gives effective, safe and painless skin tightening.



Effect

Strengthen and tighten relaxed skin with remarkable tightening effect


Effectively shrink pore, make skin texture more delicate and glossy

Remove the traces of age such as dullness and roughness, and let skin rejuvenate

Only one time care is needed to smooth wrinkles and improve them comprehensively.

Filling in shrinkage and loss of collagen vacancies to enhance the relaxation profile

Make skin tender, white and youthful



Application Scope


Type	High Intensity Focus Ultrasound
Frequency	4MHz
Power	200W
Treatment Depth	DL1.5mm DL3.0mm DL4.5mm (Standard) DL6.0mm DL8.0mm DL10.0mm DL13.0mm DL16.0mm (Optional)
Treatment Length	5 -- 25mm
Width	0 -- 10mm
Energy	0.2 -- 2.0j
Pitch	1.0mm-5.0mm
Clearanceh	1.0mm-5.0mm
Power Source	AC 110V~220V~50Hz/60Hz
Machine Size & Net weight	400X350X335mm & 8.5Kg
Package Size & Gross weight	525X445X465mm & 12Kg
Trolley Size & Net weight	540X480X790mm & 10Kg
Trolley package Size & Gross weight	680X680X330mm & 15Kg



System introduction



Aiming light



Main console



Handle connector



Handpiece



Handle support



Advantages

- Completely non-invasive and non-ablative to the epidermis layer
- No down time, and make-up can be wearied only in 15 min after the treatment
- Fast operation time , full face and neck treatment only takes 40-60 minutes
- Long lasting result, 2-3 year can last only with one treatment
- No pain, painless , safe and no side effects, anesthesia procedure is unnecessary before the treatment
- Flexible treatment, one single treatment achieves instant result
- High quality, low consumables



Treatment Tip Introduction



Item		Description
1	Top label	Tip transducer and relevant information
2	Treatment indication	Treatment range and treatment sensor

Treatment Tip Introduction



Type: ULTRALIFT-S



1.5mm



3.0mm



4.5mm



6.0mm



8.0mm



10.0mm



13.0mm



16.0mm



Treatment scope



Face

Full face
neck



Face

Full face
neck



Face

Full face
neck



Body

Back
chest
abdomen



Body

Chest
back
abdomen



Body

Back
leg
abdomen



Body

Abdomen
and legs



Body

Abdomen
and legs




Q & A

Who is suitable for Ultherapy treatment ?

It suit for 30 to 60 years old customer who want to reduce the pore size and looks younger.

when you will see the treatment result?

It suit for 30 to 60 years old customer who want to reduce the pore size and looks younger.



Q & A

why some clients have good effect of face lifting, but some clients have hollow cheek after treatment?

80% effect comes from the after treatment recovery, it's depend on patients skin property and recruitment of cell-repairing factor (such as DPHP)
(If patients age more than 30 years old, the collagen will loss more and more seriously, so before treatment, they need drink some Special repair of collagen ACMETEA so that can help regeneration of collagen fibers)

During these time, If Patients didn't support enough collagen timely, due to the deficiency of regeneration collagen will result in hollow cheek.



If a patient undergo Ultherapy treatment too frequently, what risk he will take ? what is the suitable interval time for Ultherapy treatment?







Usually the treatment results of Ultherapy treatment can last 2-3 years. Ultherapy treatment has a little damage on Subcutaneous, so we don't suggest a patient accept Ultherapy treatment too frequently, we suggest the interval time is above 2 years.

How to do the postoperative care after Ultherapy treatment?

After the treatment, the patient can do the daily life. There is no special considerations. The treatment area maybe swelling, but don't worry, it will disappear after a few hours.



Operator interface



Treatment procedure

Step 1



Deeply clean the skin

Step 3



Draw the line confirm treatment area

Step 2




Smear Ultrasound Gel

Step 4



Choose good parameter doing treatment accordingly

Before and After Treatment



THANKS

Welcome to consult more!

www.sahemedical.com

Address: 2nd floor, Building No.7, Maohuagongchang, Nancaizhen, Shunyi District, Beijing

Tel: 86-15028586078

Email: emma@sanhebeauty.net

